

# LABORHITS

## Dry Ice Machine SnowPack®

**Ideal coolant preparation for sample transport**

**Easy to handel**



### Feature

- Produces dry ice tablets at  $-79\text{ }^{\circ}\text{C}$  with 50 or 250 g
- Operates without electricity
- Taint and odour free



The dry ice machine SnowPack® produces dry ice tablets of  $-79\text{ }^{\circ}\text{C}$  with 50 or 250 g weight within a short time. Suitable for all European siphon carbon dioxide cylinders with valve without pressure regulator (W 21.80 x 1¼" acc. to DIN 477, no. 6).

- Gassing period (per tablet): approx. or 1 min
- Tablets per 22 kg bottle (dependent on ambient temperature and pressure): either 90 or 16 pieces approx.

**Applications:** hospitals, universities, schools, laboratories, chemical industry, materials testing.

**Examples of use:** lowering the temperature in cooling baths or cold traps, pre-cooling tissue samples, ensuring safe transport or embedding biological substances, freezing blood and tissue samples, studies in chemistry, physics and materials technology, etc.

Diam. mm	Height mm	Tablet weight g	PK	CatNo.
50	30	50	1	<b>9.253 006</b>
80	60	250	1	<b>9.253 008</b>