

# EBULLIOMETER model Malligand ETH

## Quick method for alcohol detection



## Perfect combination of TRADITION and INNOVATION

### TRADITION

The frame keeps the original technical and productive features, as a guarantee of a perfect thermodynamic balance

### ALCOHOL DETECTION %

6 minutes

### POWER

52 working hours

Rechargeable battery with micro USB

### TECHNICAL FEATURES

- Measuring range: 5 - 20 V/V
- Unit of measure: 0,01% V/V  
(Measures in ebulliometric range)
  
- Repeatability: + 0,10
- Tolerance: + 0,10



### INNOVATION

Electronics detect the boiling temperature, calculate the thermal balance point and show the alcoholic value on the display

**EXACTA+OPTECH LABCENTER S.p.A.**

Via Bosco, 21 - 41030 S. Prospero (MO) - Italy

eno@exactaoptech.com

www.exactaoptech.com